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FTS-B265-6  
Issue 6 : Aug. 2008

## S325 Precision Cleaver Operating Instructions

- Please read entire manual prior to usage.
- This manual must be kept with the S325 Precision Cleaver.

### Safety Instructions

This manual contains the complete operating and maintenance instructions for the S325 Precision Cleaver. Please review this manual carefully before operating the S325 Precision Cleaver.

The following safety instructions must be observed whenever this product is operated, serviced, or repaired. Failure to comply with any of these instructions or with any precaution or warning contained in this manual is in direct violation of the standards of the design, the manufacturer, and the intended use of the instrument. THE FURUKAWA ELECTRIC CO., LTD. assumes no liability for the customer's failure to comply with these requirements.

### Safety Instructions

#### CAUTION

Failure to comply with any of the instructions, which are indicated by this symbol, may cause serious injury or damage to the machine.

#### ATTENTION

Failure to comply with any of the instructions, which are indicated by this symbol, may cause the device to function improperly and degrade the performance.

Please contact THE FURUKAWA ELECTRIC CO., LTD. with any questions relating to the description of this manual.

In no case will THE FURUKAWA ELECTRIC CO., LTD. be liable to the buyer, for any consequential or indirect damage which is caused by product failure, malfunction, or any other problem.

#### CAUTION

- Avoid soaking the S325 with water. Doing so may cause malfunction.
- Do not drop, throw, and handle it roughly. Doing so may cause machine's malfunction.
- Do not use in a humid, dirty, or unstable site. Doing so may cause fire, electrical shock, or malfunction.
- Please take care when disposing of fibers from the fiber waste bin. Scattered fibers may cause injury.
- The fiber waste bin for loose fiber must be installed into the S325 during operation, otherwise there is danger for loose fiber to scatter.

#### ATTENTION


- Use an appropriate tool when loosening or tightening the set screw to fix the single fiber adapter to the unit and gently tighten the screw.
- Strip the fiber so that 14 to 25mm is exposed passed your aimed cleave length. This will allow the automatic fiber waste disposal feature to work properly. For instance :  
In case of 16mm cleave length, remove the coating more than 30mm(=16+14), and less than 41mm(=16+25).
- Please clean the fiber with alcohol after stripping the coating material that was stripped by alcohol before setting it on to S325.
- Do not use Organic solvent (e.g. Acetone) to clean the rubber pads.

# FILEL

## 1. Components

Part	Part No.	Quantity
Main body	S325A-01 or S325S80-01	1
Soft case	S325A-02	1
Waste bin	S325A-03	2 (1 Reserve)
Single fiber adapter	S325A-04	1
Big waste box	S325A-05	1
Operating Instructions	FTS-B265	1

## 2. External Description

## 3. Applicable fiber type for S325A


Type	Cleave length	Used holder or adapter
Single fiber Coating diameter 0.25 or 0.9mm	5 to 20mm	Supplied single fiber adapter, S325A-04
Ribbon fiber 0.25mm pitch 2 to 12 fiber	10mm	Optional ribbon fiber holder

S325S80 is a cutting instrument cleaver designed only for single fiber whose optical fiber diameter is 80  $\mu\text{m}$ . Please use S325A for regular 125  $\mu\text{m}$  single fiber except this and cutting of ribbon fiber cleaving task.

## 4. Operating instructions

### 4.1. Cleaving fiber

- 1) Open lid by pressing lid release lever. Place a prepared fiber (stripped and cleaned) across cleaving area with fiber end inserted into the waste bin. Ensure that the stripped fiber is long enough to reach the pinch roller.  
In the case of ribbon fiber cleaving, please use the ribbon fiber holder (Optional parts).  
When using ribbon fiber, check that the fibers are aligned and no fibers overlap with each other.
- 2) Close the lid and push the lid release lever until it stops. By pushing it completely, the blade will slide to cleave the fiber and send the loose fiber to the waste bin automatically.  
Then the lid will automatically release.
- 3)  Slide the sliding blade slowly to prolong the blade life.
- 4) Open the lid and take out the fiber from the main body.

## FILE

### 4.2. Disposing of loose fiber


Dispose of loose fiber as follows after completion of work and before transporting.

- 1) Grab the sides of the waste box with fingers and pull it.
- 2) Grasp the upper part of waste bin with fingers and open the lid while squeezing lid lightly.  
Turn waste box upside down and dispose of loose fiber in a trash can.

Do not put finger inside of waste bin because pieces of the optical fiber sometimes will remain in waste box after emptying.

Return waste box to original position after disposing of loose fiber.

**⚠ATTENTION** : When inserting waste box into main body, press waste box until click is heard.

### 5. Maintenance and Rotation of blade

The blade has 16 positions printed on the side of the blade (i.e. "0-2-...-14-"). Each position can cleave up to 2,000 cleaves. However, cleaving amount may be decreased by dust or oil. When cleaving quality degrades, clean the top of blade and upper and lower bending pad using a cotton swab soaked with ethyl alcohol.

**⚠ATTENTION** : Do not use Organic solvent on rubber (e.g. Acetone).



**⚠ATTENTION** : Please note 2,000 cleaves per position is not guaranteed. The number of cleaving changes depending on the situation and environment (e.g. dust and oil.)

If above operation does not improve the cleave quality the blade needs to be rotated. Follow the below procedures to rotate the blade.

- 1) Remove the waste bin from the main body. Loosen the blade set screw using a screwdriver.
- 2) Use a cotton swab, etc. to rotate the blade so that the number shown at the top of the blade increase by one step. (Counter clock wise)  
The blade is positioned at "0" when delivered from factory. Next position is "-" right of "0".

**⚠ATTENTION** : Do not touch the blade directly. It is very sharp.

- 3) Tighten the blade set screw again.

After all 16 positions are exhausted, the blade must be replaced by new one at customer's expense. Contact service center for replacing blade.

## FITEL

This cleaver has an automatic waste fiber disposal function. To keep it operating properly, clean the pinch roller with cotton swab soaked in alcohol after use.


After repeated use, abrasion may form on the rubber of the pinch roller. However, the pinch roller keeps the function of collecting optical fiber, as long as the depth of a groove is smaller than the diameter of the fiber. When the function cannot be recovered, even if after cleaning the pinch roller with alcohol, replace the pinch roller with a new one at customer's expense.

If the ring ring does not move smoothly, spray commercial lubricant near the rotating shaft connecting the blade to the release lever. Take care not to spray the lubricant in other areas. Then push the lid release lever several times to apply the lubricant thoroughly.

### 6. Operating with ribbon fiber holder

Please go after detaching remove single fiber adaptor (S325A-04) from the main body when you use the ribbon fiber holder (Optional parts).

### 7. Contacts

## THE FURUKAWA ELECTRIC CO., LTD.

#### Head Office :

2-3 Marunouchi 2-chome, Chiyoda-ku, Tokyo, 100-8322 Japan  
Interconnectivity Equipment, FITEL Products Division, Telecommunications Company  
TEL : +81-3-3286-3445 FAX : +81-3-3286-3708

#### Sales & Marketing Global Business Solutions, Telecommunications Company

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E-MAIL [formulary@montana.com](mailto:formulary@montana.com)

### JD-4 (JARB) HOLOGRAPHY PROCESSING KIT

**All the chemicals needed for making holograms using Slavich PFG-03M plates.**

Holographic film plates, supplies, and instructions for holography are available from:

**Integraf LLC**

Holography Supplies & Resources

(425) 821-0772

[Http://www.integraf.com](http://www.integraf.com)

The JD-4 holography developer and processing kit provides all the chemicals needed for making holograms using Slavich PFG-03M or ColourHolographics BB-640 holographic emulsions. Simply mix the dry chemicals in the JD-4 kit with water to prepare the developer components and bleach solution. This formula has been proven to be most useful for all kinds of holograms. It is relatively safe, even for the home hobbyist.

The optimization of using JD-4 for quick processing of holograms recorded on Slavich PFG-03M plates and film was made by Tung H. Jeong, Riley Aumiller, Raymond Ro, and Jeff Blyth; thus it is called the **JARB** processing regime.

JD-4 is a convenient alternative for the GP-2, a developer that has been used for years for making reflection holograms on PFG-03M emulsion. GP-2 has certain limitations, including its long processing time of 12-15 minutes and its natural drying time of over an hour. GP-2 also does not address the fact that the PFG-03M emulsion requires long exposure times (ten times more than PFG-01, for example), and is extremely soft, like gelatin, and often floats off the glass during development and washing.

The introduction of JD-4 increases the sensitivity of the PFG-03M 1.5 mJ/cm<sup>2</sup> (milli-joules per square centimeter) to 0.15 mJ/cm<sup>2</sup>. This therefore shortens the exposure time for holograms to one-tenth as long as when processed in GP-2. JD-4 also hardens the emulsion so that the holographic plates is easier to develop and can be dried using warm air, such as that from a hair dryer. As a result, the total processing time using JD-4 can be as short as three minutes, from developing to drying!

Given the above, JD-4 is ideal for making holograms during a lecture demonstration or for laboratory exercises or workshops where many students must make holograms in a limited time.

### FOR YOUR CHEMICAL SAFETY

Like many household cleaners, chemicals in general should be considered dangerous and must be treated with respect. Please read all the warning labels on each package. It is good practice to use eye goggles, dust mask, apron and rubber gloves when mixing chemicals. While the chemicals have low volatility, working in a ventilated area is recommended. Although most chemicals in JD-4 are considered non-hazardous by the EPA, the kit does contain small amounts of chemicals that the EPA does consider hazardous.

**SODIUM HYDROXIDE:** Although your kit contains only a small amount of solid Sodium Hydroxide, it must be treated with special care. Sodium Hydroxide, as a solid or in solution, is a dangerous chemical. It is a corrosive and when spilled on the skin, will cause a chemical burn. Its action is insidious because the burn occurs without pain. When working with Sodium Hydroxide wash your hands frequently and without soap. If you detect a soapy feeling while washing, Sodium Hydroxide is present. In such a case, wash thoroughly with soap and water. Beads or pellets of solid Sodium Hydroxide are easily spilled during solution preparation. If spillage occurs outside of a sink, all of the spilled solid must be cleaned up, use a damp disposable towel. If the solid is not cleaned up, it will absorb the moisture from the air and form a puddle of very caustic hydroxide, which will not evaporate. We strongly urge you to wear both safety glasses and rubber gloves when working with solid Sodium Hydroxide or its solutions. We strongly urge you to wear both safety glasses and rubber gloves when working with solid sodium hydroxide and its solutions.

**METOL:** Some individuals become sensitized (develop allergic symptoms or rashes) when using metol. If this should occur, discontinue use and consult a physician.

All other chemicals contained in this kit are considered non-hazardous, but we ask you to still use care by using a dust mask and rubber gloves.

The user assumes all risks upon accepting these chemicals. IF FOR ANY REASON YOU DO NOT WANT TO ASSUME ALL RISKS PLEASE RETURN THE KIT WITHIN THIRTY (30) DAYS FOR A FULL REFUND.

Consult with local sewer and water authorities regarding proper disposal of darkroom chemicals in your area.

## MIXING THE STOCK SOLUTIONS

The chemicals in this JD-4 kit are already pre-measured to make approximately 1000ml of each of three solutions. You will need to add de-ionized or distilled water, which can be purchased at your local grocery store. It is best to use distilled water that contains no other chemicals although distilled drinking water, which sometimes contains small but negligible amounts of other chemicals, can also be used. Water from your tap generally contains fluoride and other impurities that may reduce the quality of your hologram. For efficiency and safety, teachers or an adult should pre-mix the solutions for students.

### PART A SOLUTION

CHEMICAL	AMOUNT	AMOUNT
	(Cat. no. 04-3040)	(Cat. no. 04-3045)
Metol or Elon (p-Methylaminophenol sulfate)	4 grams	40 grams
Ascorbic Acid (powder)	25 grams	250 grams
<b>ADD</b> Distilled water	1000 ml	10000 ml

### PART B SOLUTION

CHEMICAL	AMOUNT	AMOUNT
Sodium Carbonate, Anhydrous	70 grams	700 grams
Sodium Hydroxide	15 grams	150 grams
<b>ADD</b> Distilled water	1000 ml	10000 ml

### BLEACH SOLUTION

CHEMICAL	AMOUNT	AMOUNT
Copper Sulfate, Pentahydrate	35 grams	350 grams
Potassium Bromide	100 grams	1000 grams
Sodium Bisulfate, Monohydrate	5 grams	50 grams
<b>ADD</b> Distilled water	1000 ml	10000 ml

Use three 1 liter (or larger) size clean glass or plastic bottles with leak-proof caps. Label them **Part A**, **Part B**, and **Bleach** respectively. To help dissolve the chemicals, you can heat the water until it's lukewarm. But don't use the solution until it returns to room temperature. Optionally, you can also prepare each solution in a clean beaker and then pour the solution into the bottle.

**PART A.** Fill the bottle marked Part A with roughly 1000 ml of distilled water. Dissolve the metol first. Then add the ascorbic acid. Tightly cap the bottle. Part A will oxidize if it is exposed to oxygen. In time (over a few weeks to a few months), the solution may turn yellow due to the oxidation of the ascorbic acid – the solution is still usable. Once the solution turns dark brown, the potency is lost and it must be discarded.

#### Helpful Tip

To extend the shelf-life of Part A, protect it from oxidation. To do so, subdivide the solution into smaller bottles so that the unused portions are in full-capped bottles, with little or no air space. For example, try using two 500ml bottles for Part A instead of one 1000ml bottle. Refrigeration also slows down oxidation. Do not refrigerate Part A with food to prevent possible mistaken identity as food. When stored appropriately, the solution can last for months.

**PART B.** Follow a similar procedure as above for Part A, dissolving each chemical one-by-one, in any order. This solution will keep for many months at room temperature.

**BLEACH.** Follow a similar procedure as above for the Bleach, dissolving each chemical one-by-one, in any order. This solution will keep for many months at room temperature.

Store all chemicals in a safe place away from food and children.

## HOLOGRAM EXPOSURE

For detailed instructions on making holograms, study the article "Simple Holography" found on Integrat's website ([www.integrat.com](http://www.integrat.com)). Before making your exposures, you should mix and prepare your chemicals as follows.

## PREPARATION

Have the following items on hand:

- Your pre-mixed JD-4 Part A, Part B and Bleach Solution.
- 1 additional gallon (4 liters) of distilled or de-ionized water for best results. Tap water will also work here but not as well. Avoid hard water. This will be used to rinse the holograms between each processing step.
- 2 small glass or plastic trays, just large enough so that the hologram you are making can be submerged in a horizontal position.
- 2 large glass or plastic trays to hold 1 liter of distilled water for rinsing. Tap water may also work but not as well.
- 1 (optional, but recommended) large tray to hold 1 liter of distilled water mixed with about 1 ml of photographic wetting agent such as PhotoFlo or Formaflo. You can also use a small tray with less wetting agent, but you should replace the solution after a few holograms.
- 1 rubber glove

Now label one small tray as Developer A&B. Then, mix equal portions of Part A and Part B, enough so that the hologram to be developed can be totally submerged. Once mixed, the combined A&B solution can be used to develop several holograms, and can last several hours.

Next to the developer tray place a large tray with one liter of distilled water. This will be used as a rinse.

Next, label another small tray as Bleach. Put enough bleach into it so that the hologram can be totally submerged.

Next to the bleach place another large tray with one liter of distilled water. This will be used as a rinse.

Optionally, place a large tray with the wetting solution in 1 liter of distilled water. Using a wetting solution is optional but recommended. It allows the hologram to dry evenly, thus helping you prevent smudges or streaks.

Check the order of the trays: developer A&B, rinse, bleach, rinse, wetting solution.

## PROCESSING PROCEDURES

After the holographic plate is exposed, hold it by the edges with your glove hand (or tongs). Keep the emulsion (sticky) side facing upwards to protect the emulsion from accidentally scraping the bottom of your developing tray. Complete the following steps in a dark room. You can use a green safelight. Alternatively you can use a standard night light without allowing any direct light to the holographic plate. (After the bleaching process, it is safe to turn on the lights, if preferred.)

### 1. Develop:

Quickly submerge the plate into the developer so that all parts get wet evenly. Slush it around for about 20 seconds until the holographic plate turns uniformly black.

### 2. Rinse:

Rinse the developed hologram with slight agitation for at least 20-30 seconds. For best results and longer lasting holograms, rinse up to 3 minutes to make help ensure that all of Part A has been rinsed off.

### 3. Bleach:

Place the rinsed hologram into the bleaching solution; agitate it until the plate is completely clear (this may take up to 1 minute); bleach for another 10-20 seconds.

### 4. Rinse again:

Rinse the bleached hologram with agitation for at least 20-30 seconds (up to 3 minutes).

### 5. Finish in wetting solution:

Optionally, place the finished hologram in this solution for about 20 to 60 seconds. Then, remove the hologram to dry. For best results, avoid streaks or runs as you remove the hologram from the solution. You can actually do this step under light so you observe if the wetting solution is evenly coated. (After the bleaching process, the hologram is safe to process under regular indoor lighting).

A good way to dry the hologram is to stand it against a vertical surface with the bottom edge resting on a hand-towel or tissue paper. Best results are obtained when it dries naturally in clean, dust-free air. However, if time is limited, the hologram can be quick-dried by holding it vertically and blowing warm air across it with a hair dryer.

For a reflection hologram, the image can be viewed after thorough drying, which may take minutes to hours, depending on ambient conditions and technique. Transmission holograms, on the other hand, can be viewed with a laser even when wet.

**Keep from children. Always label packages clearly as poisonous materials not to be consumed.**

## PHOTOGRAPHERS' FORMULARY INC.

PO Box 950 • Condon MT 59826 • 406-754-2891 • FAX 406-754-2896  
E-MAIL [formulary@montana.com](mailto:formulary@montana.com)

### PHOTOGRAPHERS' FORMULARY SAFETY BULLETIN: HYDROXIDES

Hydroxides that are commonly used in photographic practice include Ammonium Hydroxide solution ( $\text{NH}_4\text{OH}_2\text{O}$ ), Lithium hydroxide ( $\text{LiOH}$ ), Potassium Hydroxide ( $\text{KOH}$ ), and Sodium Hydroxide ( $\text{NaOH}$ ). Like many common household cleaners, your JD-4 kit contains sodium hydroxide, but in a very small amount (only 15g).

Each of these chemicals, in solid or liquid form, is extremely caustic. Since caustic materials by definition are capable of dissolving protein, including animal tissue, one should understand the behavior of these materials and the proper techniques for handling them.

Solutions of the hydroxides, if spilled on the skin, will slowly dissolve it, and if splashed in the eye, can cause blindness in a short time. The dry material is hygroscopic, and will absorb water from the air or body to form a caustic liquid very readily.

**IF HYDROXIDES CONTACT THE SKIN:** Wash the area thoroughly with running water until the slipperiness is no longer present. Slipperiness is due to the hydroxide dissolving the skin. If desired, rinse with vinegar or working strength acid stop bath. Wash thoroughly. Treat damaged tissue as a burn.

**IF HYDROXIDES CONTACT THE EYE (dry or wet):** Place your head at once under running water (the sink is fine) and wash the eye for 5 to 10 min. Don't bother with eyewashes, etc. Time is all-important. **WASH THE EYE FIRST**, then call a doctor at once.

To dissolve the hydroxide, simply stir the pellets into the solvent. It will dissolve very readily. It will not be necessary to pulverize the pellets or flakes. Large amounts of heat are liberated when hydroxides are dissolved, and if care is not taken, glassware may be broken or spattering may occur. It is prudent to dissolve a portion of the hydroxide, allowing the solution to cool before proceeding. Use cold or ice water when dissolving hydroxides in general. (For JD-4, the amount of hydroxide in the kit is not sufficient to create noticeable heat and thus may not require cold water.)

Be sure to pick up all the pellets that accidentally spill. The solid material will pick up moisture from the air and in time, a very concentrated, and very caustic solution forms. Dispose of hydroxide solutions by flushing down the sink or toilet with large amounts of water.

Since damage to flesh or eye can be serious, we strongly suggest the use of safety glasses and gloves when handling caustics. The use of beakers with handles is advantageous, as they are less likely to be dropped. The most important safety precaution is to take the time to move deliberately and carefully. Caustics should not be handled in the presence of children or pets.


All substances can be dangerous. Any material can be handled with safety if the correct precautions are followed. In many years of handling caustic solutions and other potentially hazardous chemicals, we have had no serious difficulty, and with a reasonable amount of care, you need have no problems. We counsel respect, but not fear.

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## COMPANY PROFILE

Lambda Scientific Systems, Inc. specializes in developing and marketing scientific instruments and systems designed and manufactured specifically for experimental education in physics at colleges and universities. Our mission is to become a premier supplier of high-quality, robust, easy-to-use, and affordable scientific instruments and systems to college educators and students for their teaching and learning of both fundamental and modern experiments in physics. Our products focus on comprehensive physics education kits, as well as light sources and opto-mechanic components.

Our physics education kits cover a wide range of experiments in general physics, especially in geometrical optics, physical optics, and fiber optics. Experiments include lens imaging, interferometry, diffraction, holography, polarization, laser physics, quantum optics, and Fourier optics through a series of the most representative apparatus such as Newton's ring apparatus, Young's modulus apparatus, Michelson interferometer, Fabry-Perot interferometer, Twyman-Green interferometer, Fourier spectrometer, and laser.

Our fiber optics education kits keep a pace with the advent of fiber optical communication technology by designing experimental systems to teach fundamental optical fiber concepts such as fiber-to-fiber coupling, fiber-to-source coupling, fiber numerical aperture, fiber mode, fiber transmission loss, and fiber sensing. These kits also give students an opportunity to be familiar with modern fiber optic components or apparatus such as Mach-Zehnder interferometer, variable optical attenuator, fiber isolator, fiber splitter, fiber switch, wavelength-division multiplexer, fiber amplifier, and transmitter.

Our light sources include Xenon lamp, Mercury lamp, Sodium lamp, Bromine Tungsten lamp and various lasers.

We also provide a variety of opto-mechanical components such as optical mounts, optical breadboards and translation stages. Our products have been sold worldwide. Lambda Scientific Systems, Inc is committed to providing high quality, cost effective products and on-time delivery.

## Contents

1. Room Light Holography with White-Light Reconstruction.....	1
2. Fresnel (Transmissive) Hologram.....	1
3. Reflective Hologram.....	2
4. Image Plane Hologram .....	3
5. Two-Step Rainbow Hologram .....	5
6. One-Step Rainbow Hologram.....	6
7. Laser Safety .....	7
8. Part List.....	7
Appendix I User Instructions of Exposure Timer.....	9
Appendix II User Instructions of Photopolymer Plate .....	10

## 1. Room Light Holography with White-Light Reconstruction

The so-called “room light operation” refers to as holographic photography and development of photopolymer holographic plates under room light. In addition, the obtained hologram can be reconstructed with while light. As compared to silver salt plates that must be operated in a dark room for photo taking and developing, room light operation of photopolymer plates is more convenient with higher diffraction efficiency. Typically, the diffraction efficiency of silver salt plates does not exceed 30% or 40%; while the diffraction efficiency of photopolymer plates can be as high as 80% or 90%.

The holographic photography of a photopolymer plate is a photochemical process during which micro molecules or monomers polymerize into macro molecules or polymers. Thus, the obtained hologram is a phase hologram; while a silver salt plate is exposed to light, its silver halide particles undertake photochemical reaction and accumulate potential energy by absorbing optical energy to form metal silver particles. After developing, the light-exposed silver particles deposit and stay while other silver halide particles are washed away after fixing process. The obtained hologram is an amplitude hologram.

Although the sensitivity of photopolymer plates is lower than that of silver salt plates, other attractions of photopolymer plates such as room light operation, high diffraction efficiency, and easy handling make them preferred in many applications. There are many methods for enhancing the sensitivity of a photopolymer plate, such as to raise the laser power, improve the reflectivity of the object, select stable workbench, or increase the exposure quantity.

## 2. Fresnel (Transmissive) Hologram

### Schematic of Experimental Setup


Figure 1 Schematic of experimental setup 1


Figure 2 Schematic of experimental setup 2

Laser: Semiconductor laser with 2-D adjustable holder (LEPO-20) on magnetic base (LEPO-4)

K: Exposure timer (LEPO-30) on magnetic base (LEPO-4)

BS: Beam splitter with plate holder (LEPO-13) on magnetic base (LEPO-4)

M: Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

M<sub>1</sub>: Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

M<sub>2</sub>: Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

L: Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

$L_1$ : Beam expander ( $f=6.2$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

$L_2$ : Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

H: Photopolymer holographic plate with plate holder (LEPO-13A) on magnetic base (LEPO-4)

O: Small object on object loading platform (SZ-20) with magnetic base (LEPO-4)

### **Experimental Procedure**

1. Align laser, exposure shutter, beam splitter, beam expanders, mirrors, sample loading platform, plate holder, and small object to be axial (at identical height).
2. Turn the laser on, configure the experimental setup based on Figure 1, and make sure that the half portion of the light transmitted through the beam expander is illuminated on the object (slightly over-illuminated on the object). Light reflected by the object to holographic plate is called the object light. Another portion of the light from the beam expander is reflected by mirror M onto the holographic plate, called the reference light. The intersection angle between the object light and the reference light should be made between  $30^\circ$  and  $50^\circ$  with an intensity ratio adjusted within a range from 1:2 to 1:5.
3. To adjust the intensity ratio of object light to reference light, an adjustable attenuator can be inserted in the object or the reference light path. In Fig. 2, beam expander  $L_1$  or  $L_2$  can be moved back or forth to change the beam size, thus the intensity of the expanded beam.
4. To change the intersection angle between object and reference light at holographic plate, arrange the location of the object or the holographic plate. Mirror M in Figure 1 or mirror  $M_2$  in Figure 2 can also be moved to change the intersection angle between object and reference light.
5. Optimal exposure time depends on the intensity of the object light reflected from the object and the laser power. Due to the low sensitivity of photopolymer, larger exposure quantity should be used. One can cut a small strip off a holographic plate and expose it with different time duration. With a semiconductor laser of 40 mW, optimal exposure time should be around 90s to 120s.
6. Upon the completion of optical configuration, make sure all the magnetic based are secured on the optical table. Put a photopolymer plate on the plate holder with the film side towards the laser. Fasten the holographic plate on the plate holder and wait 2 minutes until the stress added the locking screw of the holder is released. Expose the holographic plate with desired exposure time duration.
7. For post-exposure processing procedure of holographic plates, read the user instructions of light sensitive photopolymer (RSP) plates.

### **3. Reflective Hologram**

The difference between reflective and transmissive hologram is judged by the holographic plate set as the boundary. If both object light “O” and reference light “R” lie on the same side of the holographic plate as shown in Figure 3 (a), then the formed hologram is transmissive; if the object and reference light are incident from opposite side of the holographic plate as seen in Figure 3 (b), then formed hologram is reflective.


Figure 3 (a) Transmissive hologram


Figure 3 (b) Reflective hologram

### **Schematic of Experimental Setup**


Figure 4 Schematic of reflective hologram recording

Laser: Semiconductor laser with 2-D adjustable holder (LEPO-20) on magnetic base (LEPO-4)

K: Exposure timer (LEPO-30) on magnetic base (LEPO-4)

L: Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

H: Holographic plate with plate holder (LEPO-13) on magnetic base (LEPO-4)

O: Small object on object loading platform (SZ-20) with magnetic base (LEPO-4)

### **Experimental Procedure**

1. Align laser beam and other optical components such as exposure shutter, beam expander, plate holder, and small object to be axial (at identical height).
2. Move beam expander back or forth, so the expanded laser beam is slightly over illuminated the small object (do not over-illuminate the object-losing optical energy; or under-illuminate the object-losing object image).
3. Place the object as closely as possible to the holographic plate to enhance the reflection from the object to the holographic plate.
4. Upon the completion of optical configuration, make sure all the magnetic based are secured on the optical table. Put a photopolymer plate on the plate holder with the film side towards the laser. Fasten the holographic plate on the plate holder and wait 2 minutes until the stress added the locking screw of the holder is released. Expose the holographic plate with desired exposure time duration (20s to 30s if a glass or ceramic object is used and the laser power is around 40 mW).
5. For post-exposure processing procedure of holographic plates, read the user instructions of light sensitive photopolymer (RSP) plates.

### **4. Image Plane Hologram**

The difference between Fresnel and image plane hologram is that the former is recorded with object light diffused directly from an object to a holographic plate; while the latter employs an imaging lens to form an image of an object near a holographic plate to interfere with a ray of reference light as shown in Figure 5.


Figure 5 Schematic of image plane hologram recording

Laser: Semiconductor laser with 2-D adjustable holder (LEPO-20) on magnetic base (LEPO-4)

K: Exposure timer (LEPO-30) on magnetic base (LEPO-4)

BS: Beam splitter with plate holder (LEPO-13) on magnetic base (LEPO-4)

$M_1$ : Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

$M_2$ : Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

$L_1$ : Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

$L_2$ : Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

$L_3$ : Imaging lens ( $f=225$  mm) with lens holder (LEPO-9) on magnetic base (LEPO-4)

H: Holographic plate with plate holder (LEPO-13) on magnetic base (LEPO-4)

O: Small object on object loading platform (SZ-20) with magnetic base (LEPO-4)

### **Experimental Procedure**

1. Align laser beam and other optical components such as exposure shutter, beam expander, mirrors, lenses, plate holder, and small object to be axial (at identical height).
2. Based on Figure 5, align object light path and let  $BSM_1 + M_1O + OH = A$ . Move  $L_1$  back and forth until object  $O$  is fully illuminated. Insert lens  $L_3$  between object  $O$  and holographic plate  $H$  and move  $L_3$  back or forth to image the object near the holographic plate (either at the front or at the back of the plate).
3. Align reference light path and let  $BSM_2 + M_2H = A$ . Insert  $L_2$  in the reference light path and move  $L_2$  or  $L_1$  back and forth until the intensity ratio of object light to reference light falls within the range of 1:2 and 1:5. The intersection angle between object beam and reference beam should be between  $30^\circ$  and  $50^\circ$ . If not, adjust the locations of mirror  $M_1$ , or  $M_2$  or object  $O$ .
6. Upon the completion of optical configuration, make sure all the magnetic based are secured on the optical table. Put a photopolymer plate on the plate holder with the film side towards the laser. Fasten the holographic plate on the plate holder and wait 2 minutes until the stress added the locking screw of the holder is released. Expose the holographic plate with desired exposure time duration (90s to 120s if a glass or ceramic object is used and the laser power is around 40 mW).

7. For post-exposure processing procedure of holographic plates, read the user instructions of light sensitive photopolymer (RSP) plates.

## 5. Two-Step Rainbow Hologram

Typically, an off-axis Fresnel hologram (or transmissive hologram) can only be reconstructed with laser light. The so-called two-step rainbow holography is to view a laser reconstructed hologram with white light. First, a hologram is recorded with laser light. Then, reconstruct the hologram with laser light so that the reconstructed real image of the original object is formed near a second holographic plate to be interfered with a ray of reference light. The second hologram can be reconstructed with white light, as shown in Figure 6.


Figure 6 Schematic of two-step rainbow hologram recording

Laser: Semiconductor laser with 2-D adjustable holder (LEPO-20) on magnetic base (LEPO-4)

K:Exposure timer (LEPO-30) on magnetic base (LEPO-4)

BS: Beam splitter with plate holder (LEPO-13) on magnetic base (LEPO-4)

M<sub>1</sub>:Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

M<sub>2</sub>:Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

M<sub>3</sub>:Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

L<sub>1</sub>: Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

L<sub>2</sub>: Beam expander ( $f=6.2$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

L<sub>3</sub>: Collimating lens ( $f=225$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

L<sub>4</sub>: Imaging lens ( $f=150$  mm) with lens holder (LEPO-9) on magnetic base (LEPO-4)

G: Slit screen with plate holder (LEPO-13) on magnetic base (LEPO-4)

H<sub>1</sub>: Fresnel hologram with plate holder (LEPO-13) on magnetic base (LEPO-4)

H<sub>2</sub>: Holographic plate with plate holder (LEPO-13B) on magnetic base (LEPO-1)

## Experimental Procedure

1. Based on Figure 6, align object light path with beam expander  $L_2$ , collimating lens  $L_3$ , and imaging lens  $L_4$  so that object light is incident on Fresnel hologram  $H_1$ . The reconstructed real image is formed on holographic plate  $H_2$  through slit  $G$ .
2. As object light is imaged through a narrow slit, the spectral bandwidth of the reconstructing light is suppressed by the slit when reconstructing the hologram thus avoiding a blurry image. The width of the slit should be around 5 to 8 mm. If the width of the slit is too wide, then the reconstructed image becomes blurry; if the width is too narrow, the intensity of the object light becomes inadequate.
3. For white light reconstruction, the direction of the slit should be orientated vertically to maintain the horizontal parallax of the object. When the hologram or the observation eye is moved horizontally, the horizontal parallax of the object will be observed; when hologram or eye is moved vertically, only the color of the reconstructed image varies.
4. Repeat steps 6 and 7 in sec. 2.

### 6. One-Step Rainbow Hologram

One-step rainbow holography removes the tedious two steps in two-step rainbow holography by introducing a slit in before or behind the imaging lens the optical path as shown in Figure 7.


Figure 7 Schematic of one-step rainbow hologram recording

Laser: Semiconductor laser with 2-D adjustable holder (LEPO-20) on magnetic base (LEPO-4)

K: Exposure timer (LEPO-30) on magnetic base (LEPO-4)

BS: Beam splitter with plate holder (LEPO-13) on magnetic base (LEPO-4)

M: Plane mirror ( $\Phi 40$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

$L_1$ : Beam expander ( $f=4.5$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

$L_2$ : Collimating lens ( $f=150$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

$L_3$ : Beam expander ( $f=6.2$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-3)

$L_4$ : Collimating lens ( $f=225$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

$L_5$ : Imaging lens ( $f=150$  mm) with 2-D adjustable holder (LEPO-8) on magnetic base (LEPO-4)

G: Slit screen with plate holder (LEPO-13) on magnetic base (LEPO-4)

O: Small object on object loading platform (SZ-20) with magnetic base (LEPO-4)

H: Holographic plate with plate holder (LEPO-13B) on magnetic base (LEPO-4)

### **Experimental Procedure**

1. Align laser beam and other optical components such as exposure shutter, beam expander, mirror, lenses, slit, plate holder, and small object to be axial (at identical height).
2. Based on Figure 7, align object light path with beam expander  $L_1$  and collimating lens  $L_2$ , so that object light is slightly over-illuminate small object  $O$ . Place adjustable slit  $G$  near object  $O$  as shown in Figure 7 and adjust the width of the slit around 5 to 8 mm. Now the object light forms a real image of the object by imaging lens  $L_5$  through slit  $G$  near holographic plate  $H$ . Measure the length of the object light path as  $BSO+OH=A$ .
3. Align reference light path by positioning the location of mirror  $M$  to let light is reflected by mirror  $M$  to holographic plate  $H$ . Make  $BSM+MH=A$ . Put beam expander  $L_3$  and collimating lens  $L_4$  in the optical path so that collimated reference light is illuminated on plate  $H$ .
4. Move  $L_1$  or  $L_3$  back or forth to change the beam size of the object or reference light so that the intensity ratio of object to reference light is in the range of 1:3 and 1:5. Alternatively, this can be done by inserting an optical attenuator in the object or reference light path.
5. Adjust the location of object, or mirror, or holographic plate so that the intersection angle between object and reference light is within  $30^\circ$  to  $40^\circ$  at holographic plate.
6. Repeat steps 6 and 7 in sec. 4.

**Note:** 1. Limited by the narrow slit, the information of object light is recorded around the region of the slit so that the reconstructed image of a rainbow hologram is no longer stereoscopic along the direction perpendicular to the slit direction. In addition, a full-size image of the object cannot be reconstructed from a partial hologram.

2. When recording a rainbow hologram, the direction of the slit is orientated vertically to maintain the horizontal parallax of the object. When the hologram or the observation eye is moved horizontally, the horizontal parallax of the object will be observed; when hologram or eye is moved vertically, only the color of the reconstructed image varies.

3. Use a ground glass screen to locate the position of the real image of the object. Adjust the object distance to the imaging lens so that a real image of the object with desired image size is obtained. If a partial image of the object is seen, then the width of the slit can be increased (do not exceed 1 cm!). Otherwise, a smaller object or an imaging lens with longer focal length should be used.

## **7. Laser Safety**

**Warning:** the semi-conductor laser in this kit is a class IIIb laser. Avoid direct eye exposure to the laser beam!

Follow the corresponding laser safety guidelines based on AS/NZS 2211.1:1997 and other lab instructions about optical components etc.

## 8. Part List

Product Code	Product Name and Description	QTY
LLL-1	Semiconductor Laser	1
LEPO-20	Semiconductor Laser Holder	1
LEPO-30	Exposure Shutter and Timer	1
	Fixed Ratio Beam Splitter (5:5 & 7:3)	1 each
	Photopolymer Holographic Plate	1 box
LEPO-13	Plate Holder	2
LEPO-13B	Plate Holder (with loading stage)	1
	Tri-Color Safety Lamp	1
	Lens (f=4.5 mm and 6.2 mm)	1 each
	Lens (f=150 mm)	2
	Lens (f=225 mm)	1
	Plane Mirror	3
LEPO-4	Universal Magnetic Base	10
	Continuously Variable Beam Splitter	1
LEPO-9	Lens Holder	2
LEPO-8	Two-Axis Adjustable Holder	6
SZ-20	Sample Stage	1
	Small Object	1
	Electric Dryer	1
SZ-43	Ground Glass	1
	Small White Screen	1
LEPO-3	Height Adjustable Holder	2
LEPO-2	Two-Axis Translation Base	1
	Illuminometer	1
	Slit Screen (rectangular hole)	1

## Appendix I User Instructions of Exposure Timer

### 1. Parts

Control box

Shutter

Joystick

### 2. Specifications

4-bit digital display

Shutter release time: 0.1-999.9 s

Time error: < 1 ms

Four modes selectable: B-gate, T-gate, Time-setting and Open

### 3. Operation

The exposure timer is composed of a photo-gate, a timer, and a joystick. Connect the power cable and photo-gate to the back panel of the control box and plug the joystick in the front panel. Switch on the power supply and the LED digital tube will light up, it shows that the exposure timer is working.

#### Working modes

“Shutter”—When it is pushed down, the photo-gate will always be open, “Time Setting” button and joystick will not function. When the button is released, the photo-gate will always be closed and the “T” or “B” mode can be selected.

“T”—When the button is pushed down, the timer will work on the “T” mode. When you push the button of the joystick on time, the photo-gate will open. Push it again the photo-gate will close. The interval time of the two actions will be shown on the LED digital tube.

“B”— When the button is pushed down, the timer will work on the “B” mode. Push it down and hold on, the photo-gate will open and the timer will begin counting time until the button is released. The display time shows the duration of the pushing of the button.

Except above three modes, the timer can work on a special mode, a so-called count-down mode. This mode is available when all the buttons are released. User can use the time setting function to set the exposure time, then push down the button of joystick and the photo-gate will be open. The timer begins counting down until the value goes back to zero.

**Note:** The photo-gate must be orientated vertically.

## Appendix II User Instructions of Photopolymer Plate

Red sensitive photopolymer plate is a phase-type holographic recording medium with high diffraction efficiency and high resolution. Its photo sensitivity is between dichromated gelatin and silver salt, and exceeds dichromated gelatin by one or two order of magnitude. When illuminated with light, red sensitive photopolymer undertakes photo-chemical reaction. During the polymerization process, micro molecules or monomers polymerize into macro molecules or polymers. Because the polymer is caused by light illumination, it is called photopolymer.

### 1. Specifications:

**Wavelength:** 632.8 nm or 647.1 nm (red sensitive)

**Diffraction Efficiency:** >80%

**Resolution:** >4000 line pairs/mm

**Thickness:** 10 ~ 20  $\mu\text{m}$

**Cross Sectional Size:** 90x240 mm

**Storage:** > 1 year under ambient condition (recommended for storage under low temperature environment such as in a refrigerator)

**Safety Lamp:** under room-light operation

### 2. Handling and Processing Procedures

#### 2.1 Cutting of an unexposed photopolymer plate

The photo sensitive side of a red sensitive photopolymer plate is little sticky and there are two plastic sticks used for separating two adjacent photopolymer plates from adhesion in the container. Take a photopolymer plate with two plastic sticks from the container and place the film side of the photopolymer plate on the two plastic sticks. Use a glass cutting knife to cut the photopolymer plate into desired sizes. After cutting, if the film side is still connected, then cut the film gently with a knife.

**Note:** make sure to cut the glass substrate side first. Never cut the film side of a photopolymer plate with a glass cutting knife.

#### 2.2 Post-exposure processing of photopolymer plates

##### **Preparation of Development Solutions:**

The normal temperature for developing the plate is at 20 °C.

Please use five beakers to prepare the following solutions and label them properly.

Beaker 1: 300 ml distilled water (0% isopropyl alcohol)

Beaker 2: 120 ml isopropyl alcohol and 180 ml distilled water (40% isopropyl alcohol)

Beaker 3: 180 ml isopropyl alcohol and 120 ml distilled water (60% isopropyl alcohol)

Beaker 4: 240 ml isopropyl alcohol and 60ml distilled water (80% isopropyl alcohol)

Beaker 5: 300 ml isopropyl alcohol (100% isopropyl alcohol)

##### **Development:**

- 1) Immerse the plate in Beaker 1 for 15-30 seconds
- 2) Immerse the plate in Beaker 2 for 60 seconds
- 3) Immerse the plate in Beaker 3 for 60 seconds

- 4) Immerse the plate in Beaker 4 for 15 seconds
- 5) Immerse the plate in Beaker 5 for approximately 3 minutes until a bright faint red or yellow green image appears.

**Note:** Gently shake the solutions during all above steps. Wear rubber gloves for preventing skin from contacting with chemicals.


- 6) Remove the plate from the beaker and dry it with hairdryer until a bright yellow image is seen on the plate.

**Note:** If the image has disappeared after a certain amount of time, it is possible to bring the image back. Just re-develop the photopolymer plate as per above steps and the image should reappear.

### **Sealing the Plate**

The recorded image on photopolymer is short shelf-life and may disappear a day later if it is unsealed or not sealed properly. To seal the film between two glass plates, follow these steps.

- 1) Peel a strip of photopolymer material along the edge of the plate, as shown in the figure at the right corner.
- 2) Put some glue on the peeled area.
- 3) Cover a prepared glass plate (1 to 3 mm thick and same size as the photopolymer plate) on the emulsion side of the photopolymer plate
- 4) Clamp the two plates tightly and wait until the glue dries.

### **2.3 Troubleshooting**

<b>If You Notice</b>	<b>Possible Cause</b>	<b>You Should</b>
Dark final image	1. Under/over exposure 2. Vibration influence 3. Over blow by hairdryer	1. Increase/decrease exposure time 2. Get rid of vibration sources 3. Stop drying immediately when you see a bright yellow image
No or partial final image	Object unlighted or not illuminated completely	Check if the laser completely illuminates the object
Photopolymer comes off the substrate	Immerged in distilled water too long	Immerge for a shorter time
The image turns from red to yellow, then blue, or disappear	Over dried	Stop drying immediately when you see a bright yellow image
Ununiformed color of final image	Uneven dried	Dry the plate evenly
Image has a narrow depth of field	Low intensity of object beam or low light coherence	Place object closer to the plate or increase the coherence length of laser

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Jonathan F. Reichert, Ph.D.

### I. INTRODUCTION

"Impelled by a belief in the unity of the forces in nature, Michael Faraday sought, and in 1845 provided, the first phenomenological evidence for the connection between light and magnetism when he discovered the effect that still bears his name. He found that plane-polarized light, propagating through matter parallel to a static magnetic field, underwent a systematic rotation of its plane of polarization."

Thus begins an interesting paper by David A. Van Baak on the topic of Faraday rotation in atomic gas.<sup>1</sup> Although carefully studied from the time it was discovered, this remarkable effect was not theoretically understood until the development of quantum mechanics. References in Van Baak's paper provide some insight into current efforts both experimentally and theoretically, even for nonlinear optics. Many textbooks also discuss this topic.<sup>2</sup>

Although the effect is typically small, Faraday rotation has found important practical applications. In modern fiber optic communications networks, it is used as the canonical optical isolator.<sup>3</sup> This is an extremely important component of such networks since it prevents the multiple propagation of reflected signals. Faraday rotation has important applications in optical modulators as well as uses in analytic apparatus for mixtures of hydrocarbons.

We strongly urge the instructor to make use of these references. Copies of two American Journal of Physics articles have been included in this manual to facilitate student exploration of the literature on the Faraday effect. Although it is an important effect, it is frequently omitted in the standard courses in E&M and Optics in both the undergraduate and graduate physics and engineering curriculum.

TeachSpin originally thought of this apparatus in connection with our teaching Signal Processor/Lock-In Amplifier. It seemed to us that the Faraday effect was an ideal experiment to show both interesting physics and the power of lock-in detection for extracting weak signals from background noise. It was with this in mind that we designed FR1-A. However, in the direct DC mode, the unit can certainly be used without lock-in detection. It must be noted that this unit comes with only one monochromatic light source, a laser pointer with a wavelength of about 650 nm. Instructors who want the student to explore the interesting frequency dependence of the effect will have to supply their own polarized light sources at other wave lengths. The Polaroid film analyzer has a useful polarizing capability from 400 to 800 nm. Since the Faraday effect is a slowly varying function of wavelength, the intensity stability of the light source is more important than its frequency stability. The diode laser is an ideal source for these measurements.

## II. APPARATUS

FR1-A consists of four basic components: the Light Source, the Solenoid (magnetic field source), the Analyzer Polaroid and the Optical Detector. Each one will be described separately.

### A. The Light Source

A red laser pointer is used as the light source. It operates at a nominal wavelength of approximately 650 nm with a power output of about 3 mw. Although the output is approximately 60% polarized, the laser light is directed through a polarizing filter which increases its polarization to about 95%. Since this filter is not permanently bonded to the laser, it is important to adjust its rotation position to achieve maximum power output of the light source. This will occur when the polarization axis of the laser corresponds to that of the filter. This can easily be adjusted visually by inserting a white card at the entrance to the solenoid.

The laser requires a voltage regulated supply of 4 volts and 40 mA. ***The red banana plug should be connected to the positive terminal post of the supply. Reversing the power supply voltage will destroy the laser.*** Note: it is important to protect both leads to the laser from an electrostatic shock. Keep the connectors shorted to each other until they are connected to their supply. TeachSpin's audio amplifier unit AALS1-A has a voltage regulated (4V) supply output specifically designed for this laser diode light source.

The laser beam must be adjusted by the operator to traverse the solenoid along its central axis. This can be accomplished by manipulation of the four nylon thumb screws on the laser mount using a white card both at the entrance and exit of the solenoid. The sample glass rod, or other samples, should be removed from the solenoid during the alignment procedure. After alignment, the sample can be reinserted into the solenoid.

### B. Solenoid

The solenoid has the following physical specifications:

Length: 150 mm	Wire Size: #18 double insulated
Turns/layer: 140	DC Resistance: 2.6 $\Omega$
Number of Layers: 10	

The approximate calibration constant for the unit, measured at the center of the solenoid is:

$$B = (11.1 \text{ mT/A}) I$$

where I is in amperes and B is in millitesla.

The magnetic field does vary along the axis of the solenoid. Such variation may be significant for certain samples, particularly those which extend outside the coils. We recommend that the students measure the axial magnetic field profile using TeachSpin's Hall Effect Probe HE1-A, or other equivalent field measuring apparatus.

The maximum continuous current through the unit is 3 amperes. For times of the order of 30 seconds, however, 10 amperes can be used without damaging the solenoid or its supports. **Extreme care** should be taken to assure that the unit does not overheat when exceeding the 3 A specification.

---

Since Faraday rotation measurements can easily be performed in 10-15 seconds by a student familiar with the equipment, it is possible to use these larger solenoid currents and thus obtain larger rotations.

#### C. Analyzer Polaroid

The unit is equipped with a rotatable Polaroid film in a calibrated mount. The decal is marked in  $5^\circ$  increments. This limits the accuracy of an angular measurement to about  $2^\circ$ . This is a conservative estimate. Some students may claim  $1^\circ$  set-ability. You be the judge.

#### D. Detector

The detector is simply a photodiode connected in series to one of three resistors; 10K, 3K, 1K. The photodiode is a *current source* and is a linear photonic detector, as long as the voltage across it is less than about 0.3 volt. Saturation begins to occur when this bias voltage appears across the diode. This makes the detector nonlinear.

*It is wise to work at bias voltages well below the 0.3 volt value.* This can be accomplished by using the proper load resistor and monitoring the bias voltage. Since the laser diode light source is high intensity, it will be necessary to use the lower values of load resistors for some Faraday rotation experiments.

#### E. Sample

FR1-A comes with one glass rod sample. Its specifications are:

Diameter: 5 mm	Length: 10 cm
Material: SF - 59	

The sample can be removed from the solenoid by first removing the detector, analyzer and light source. This gives the experimenter unrestricted access to the solenoid's bore. Care should be taken not to damage the end of the glass sample. They are highly polished. Any push-rod should have a soft end covered by optical tissue or its equivalent.

### III: EXPERIMENTS

#### A. Malus's Law

With the sample removed from the solenoid, students can measure the intensity of the transmitted light through two polarizers (one on the laser, the other the analyzer) as a function of the relative angle between them. They should also derive a simple theoretical calculation of this "law."

Be Careful!! This is a simple experiment, but one must not saturate the diode detector. Use the correct load resistor.

#### B. Measure the Verdet Constant, $v$ , DC.<sup>4</sup>

The Verdet constant is defined by the relationship:

$$\theta = vB\ell$$

where  $B$  is the strength of the magnetic field in millitesla,  $\ell$  is the length in centimeters and  $\theta$ , in radians, is the angle of rotation of the axis of the polarization after passing through the medium in the magnetic field. The units of  $v$  are then: radians/millitesla centimeter. These are by no means universal units. The units given in the literature vary and include degrees/gauss centimeter.

First, you may wish to independently measure the magnetic field inside the solenoid rather than using the  $B = 11.11$  calibration we have supplied. The field inside is not constant along the entire length. It drops off to about half the maximum value at the ends. TeachSpin's HE1-A Hall Effect Probe is an ideal instrument to measure this axial field, but other Hall sensor probes will also do the job. Such a calibration will allow the student to make corrections to their determination of the Verdet constant.

Place the 10 cm glass rod inside the solenoid, supporting it on the foam rubber sleeves. Make sure it is located in the center. Using any white surface (a business card will do) align the laser beam so that the light enters and exits along the center axis of the solenoid. It is important that this be done carefully so that the entire laser beam passes through the glass rod.

Now there are several ways the Faraday rotation can be measured. Both of the methods given below should be tried by the students so that they can explain which one is more accurate and why.

##### 1. Method 1 – Extinction

- Set the analyzer polarizer at  $90^\circ$  with respect to the entrance polarizer. This gives maximum extinction of the light to the detector.
- Connect a current regulated supply to the solenoid.
- Measure the angle of rotation of the analyzer necessary to return to maximum extinction when the solenoid's magnetic field is turned on.
- Find the rotation angle for a variety of magnetic fields.

Describe what happens when you reverse the direction of the magnetic field?  
Does the magnetic field affect the photodiode detector itself?

## 2. Method 2 – Returning to Fixed Intensity Value

- a. Set the analyzer at  $45^\circ$  with respect to the entrance polarizer.
- b. Repeat the procedure described in Method 1. This time, however, adjust the analyzer angle so that the detector output returns to the value found before the magnetic field is turned on.
- c. How does the value of the initial angle affect the results?

It should be straightforward to determine the Verdet constant from either of these measurements. Which method is best? Why?

## C. Measure the Verdet Constant $AC^5$

A more accurate and more sensitive method of determining Verdet constants uses lock-in detection and ac magnetic fields. TeachSpin originally developed FR1-A to accompany its teaching Signal Processor/Lock-In Amplifiers (SPLIA1-A) since Faraday rotation was a good example of a “small” physical effect that can easily be measured with lock-in detection techniques. This was also independently realized by Jain etc. all who have written an excellent article on how this measurement is carried out. We have enclosed a copy of their American Journal of Physics article for your convenience.

Lock-In detection affords some additional possibilities for the student.

1. The relative angle at which to set the polarizer to obtain the maximum signal is easily determined.
2. The relative phase angle of the ac signal as a function of Polaroid angle can be examined.
3. The student develops an understanding of the process of obtaining an absolute measurement of  $v$  from these ac measurements.

## D. Measurements of Verdet Constants in Liquids

The ac lock-in detection method allows one to detect extremely small Faraday rotations and thus small Verdet constants without using extremely large magnetic fields. Various liquids exhibit the Faraday effect and their constants can be measured with this apparatus. A glass cell with flat optical ends and a filling tube at one end will hold an appropriate liquid sample. With such a sample holder, however, the liquid will not all be in the same magnetic field. Part of the sample will even be outside the solenoid.

This presents an interesting challenge for the student's determination of the Verdet constant. It requires that the axial field profile be measured first, and a simple computer program or hand held calculator algorithm be developed to extract  $v$  from the measurements. This can easily be done since the Faraday effect is a linear function of the magnetic field.

#### IV. REFERENCES

1. D.A. Van Baak, "Resonant Faraday rotation as a probe of atomic dispersion." *Am. J. Phys.* 64 (6), June, 1996.

---

This paper has an excellent set of references for both the historical story of this discovery and modern practical applications of the effect. The paper also contains references to student experiments and theoretical papers analyzing this phenomenon.

2. Alan Portis, "Electromagnetic Fields," Wiley, pp 467 – 474.  
Eugene Hecht, "Optics," Addison-Wesley (Third Edition), p 862.
3. Electro-Optics Technologies, Inc. [www.eotech.com](http://www.eotech.com)
4. Frank J. Loeffler, "A Faraday rotation experiment for undergraduate physics laboratory." *Am. J. Phys.* 51 (7) July, 1983.
5. Aloke Jain et. al., "A simple experiment for determining Verdet constants using alternating current magnetic fields." *Am. J. Phys.* 67 (8) August, 1999.

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